


Optimum (GL) 150 / 200 / 200M

Mode d'emploi • Bedienungsanleitung • Operating instructions

MA 185 & MA 186

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Before powering up

Danger of electrical shock !

- Do not open any electrical boxes.
- Never use a 3-wire/2-wire adapter.
- The cable sheaths must be free of defects.

Danger of pinching on the Optimum 200M !

- Do not put your hands between the back of the table and the column when the machine is on.
- Do not put your hands between the optics and the table when the power is on.

Caution


General safety precautions must be followed when operating this electrical equipment. Not following these precautions could damage the equipment or injure personnel. Safety rules vary between companies. If a conflict exists between the content of this guide and the rules of the company using this system, the more stringent rules should be applied.


Start-up

- | | |
|-----------|--|
| Base | A) Check the levelling of the machine by adjusting the levelling feet. |
| X-Y table | B) Remove the red safety plates on the sides. 2 plates for Optimum 150 and 4 plates for Optimum 200 and 200M.
C) Fit the X and Y handles on their drive axes. |
| Electric | D) Screw the scale cables to the MCA-190 and connect the MCA-190 to the computer. Connect the 110/220V plug to the mains. Turn on the switch if there is one. |
| Base 200M | E) Remove the carrying handles for the 200M. |


Optimum 150


Optimum 200

Optimum 200M


Description Optimum 150 / 200 with manual XYZ axis drives


- 1) Camera cable to the computer
- 2) Camera
- 3) Lighting control panel
- 4) Optics
- 5) Zoom position knob
- 6) Incident LED ring light (light from below)
- 7) X-axis crank
- 8) Clutch disengagement X-axis
- 9) Transmitted LED illuminator (lighting from below)
- 10) Z-axis hand-wheel
- 11) Precision column
- 12) Coaxial LED illuminator (lighting passing through the optical axis)
- 13) Table glass
- 14) X-Y table
- 15) Granite base
- 16) Machine identification
- 17) Levelling foot
- 18) Clutch disengagement Y-axis
- 19) Y-axis hand-wheel
- 20) XY or XYZ ruler box


- 21) Jack connector for scale power supply (XYZ version only)
- 22) USB2 B connector for the computer
- 23) SUB-D connector for table X-axis scale
- 24) SUB-D connector for table Y-axis scale
- 25) SUB-D connector for column Z-axis scale (optional)

Description Optimum 200M with motorized XYZ axis drives

- 26) Control display
- 27) Control panel
- 28) Main switch ON/OFF
- 29) Emergency stop
- 30) Trackball for XY axis
- 31) Keyboard
- 32) Z-axis wheel, sharpness adjustment.
- 33) Changing sectors for incident illuminator
- 34) Intensity control for incident illuminator
- 35) Intensity adjustment for coaxial illuminator
- 36) Intensity adjustment for transmitted illuminator
- 37) Light ring
- 38) Lighting switch

General use for Optimum 150 / 200 / 200M

- **Calibration** : Make sure your device is properly calibrated. To do so, please use the instructions for the software.
- **Optics 4)** : Please refer to the "[OptiVID](#)" operating instructions.
- **Software** : Use the F1 key for online help or refer to our "Quick start" or "Instructions for Use" documents.
- **Zoom 3)** : Turn the knob 5) to change the magnification. The zoom positions are indexed. You should therefore feel a "click" when you are in the required position. A sticker on the zoom allows you to know your position. Make sure you are on the same position in your software.
- **Sharpness** : Place the part to be measured on the glass 13) of the X-Y table 14). Switch to the highest magnification of your zoom 5) to get the shallowest depth of field. Use the hand-wheel 10) to adjust the sharpness of the image. You can then return to a lower magnification, the image will remain sharp.
- **Measurement** : To move your part, use the hand-wheels 7) and 19). If necessary, you can disengage the X and Y axes using levers 8) and 18).
- **Lighting 6) 9) 12)** : Adjust the light intensities using panel 3) for Optimum with manual drives. For more information, refer to the operating instructions [MA 356-041](#). In optical inspection, it is very important to control the reflections on the parts and to know the optimal conditions to guarantee the reliability of the measurement. This involves working with rather dark images. A calibration operation on a part similar to the ones you are measuring, whose exact dimensions are known to you, will make it possible to define the ideal lighting conditions.

Special use for the Optimum 200M

- **Quick start** : A quick start manual is available on our website.
- **Power On** : Press button 28). Wait for the electronics to initialize (4 seconds). Display 26) flashes. Press the C key briefly to activate the motors. The motors are now active, you can use the machine.
- **Reference socket** : C key. Be careful when activating this function, the 3 axes will start at the stop to take their references. Press the C key for more than 3 seconds to perform the reference measurement. It is strongly recommended to activate this function each time the motors are switched on.
- **Emergency stop**: If there is a risk of accident, press the red emergency stop button 29). To restart the machine, turn this button 29) to the left or right to raise it.
- **Measurement** : To move your part, use trackball 30).
- **Sharpness** : Place the part to be measured on the glass 13) of the X-Y table 14) and use the roller 32) to adjust the sharpness of the image.
- **Lights 6) 9) 12)** : Adjust the light intensities using the panel 27). The illuminators can be switched on and off independently of each other by means of switches 38). The light intensity can be adjusted by turning the knobs 34) 35) 36). Pressing these buttons activates a coarse or fine setting. The light ring 37) indicates the mode : intense = coarse / weak = fine. To change the sectors of the incident illuminator 6), use button 33). For more information on sectors, refer to operating instructions [MA 356-041](#).
- **Keyboard 31)** : See below.

Choice of travel speed.

Deactivation or activation of trackball control for X and Y axis.

Press for more than 3 seconds to save a position.
Press briefly to move the axis to the stored position.

Stop of the 3 axis X, Y and Z. To engage the axis again, briefly press the C key.

Position "0". Move the 3 axis to the reference position "0".

Taking of reference.

Troubleshooting

Descriptions :

The illuminators 6), 9) or 12) are defective.

The motor no longer works.

Camera 2) malfunctions or does not work. A "License Error" or "License Not Found" message appears.

The PC or monitor no longer works.

Table glass 13) is broken or scratched.

The crank 10) in column 11) rotates and the column isn't moving.

A clutch 8) or 18) of the table 14) no longer works.

Software problem

For any other case, please contact us by referring to the serial number on the machine identification 16)

Procedures :

1. Contact our technical service.

1. Check that the camera cable 1) is correctly inserted on the computer side.
2. Connect the camera cable 1) to another port (usually blue for USB3 models).
3. Check, by disconnecting and reconnecting it, that a device connected to the computer (external hard disk, loudspeaker, etc.) does not interfere with the camera signal.
4. Refer to our "[OptiVID](#)" operating instructions
5. If this is not enough; contact our technical department.

1. Check that the cables are correctly connected in the plugs.
2. If this is not enough; contact your IT department or, if necessary, our technical department.

1. Order a spare glass.
Optimum 150: VER-000078.
Optimum 200 and 200M: VER-007974.

1. Contact our technical support.

1. Order a replacement clutch :
Optimum 150: MEC-011889.
Optimum 200: Contact our technical support.

1. Press the F1 key on your keyboard and refer to the online help.
2. Turn off the screen saver and automatic log-out.
3. Turn off your PC completely, wait 10 seconds and turn it back on.
4. If this is not enough; contact your IT department or, if necessary, our technical support.

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